


FOR ACHES AND PAINS

IceWave™ provides relief for sore muscles, sprains, strains and other injuries due to intense exercise.

- The IceWave Patch™ from LifeWave™ is a non-invasive nanotechnology product.
- IceWave™ works immediately providing the injured area with a thermal effect.
- As shown in the actual before and after pictures to the right, the reduction of "hot spots" is clearly evident.
- Local thermal effects include the reduction of excess heat in an injured area; reduced pain, reduced soreness, and restoration of range of motion to the specific injured site.
- LifeWave™ patches have been found in studies to have excellent results with acute minor muscle spasms and injury, as well as chronic muscle spasms.



"Infrared thermal imaging studies clearly show that LifeWave IceWave Patches™ have both immediate and local thermal effects on the skin."

- Dr. Dean Clark, D.C.



IceWave™ is ideal for extreme sports including university, college, and high school athletes.

For more information on LifeWave™, please contact:

ICE - 2500

ICEWAVE
SMART TEMPERATURE TECHNOLOGY




- RELIEVES MINOR ACHES AND PAINS**
- WORKS WITHIN MINUTES**
- NOTHING ENTERS THE BODY**
- NO MESSY CREAMS OR CHEMICALS**

Software for the Human Body

When LifeWave IceWave Patches™ are placed on an individual's body over a specific region of injury, localized thermal responses occur offering an overall thermal response much like hot and cold packs. This is all done with no chemicals, creams or medication. LifeWave's patented technology is completely nontransdermal (nothing enters the body).



Infrared thermal imaging has proven that IceWave™ patches produce a thermal response over regions that are injured, strained, or in spasm. IceWave™ patches have a thermo regulating effect on the injured area similar to cold or hot packs.



Use IceWave™ patches in the place of cold or hot packs for minor sprains, strains, tendonitis, muscle soreness, aches and pain. IceWave™ can be worn during exercise and after.

Instructions for use



White Patch

goes on
the sternum




Tan Patch

goes where
the pain is



One tan patch can be placed on any place on the body where there is discomfort; however the white patch must be placed in the center of the sternum.



For Lower Back Pain

Place one white patch on the right side of the spine where the pain is, place a tan patch on the left side of the spine accordingly. Only wear one set of patches at once.

If condition does not improve within 7 days, discontinue use and see a doctor. Use product only as directed. Do not use product if you have a health condition. Should you have any concerns, please consult your healthcare professional before using.

The adhesive of this product is hypoallergenic; however, if you experience a skin irritation simply discontinue use. Should any of the contents of the patch make contact with the skin simply wash and remove. Should you experience any type of discomfort from the use of this product discontinue use. Please review all instructions and information about this product before using.
STORE PACKAGE IN A COOL DRY PLACE.

IceWave™ Gets Results

Greater than 80% of people who have been evaluated and patched have received positive results. Approximately 20-25 percent of the people received relief within a few minutes.

Testimonials:

"I rode 60 miles on my bicycle and started feeling soreness or lactic acid in the quadriceps at 28 miles, but felt recovery on the hills the last 8 miles."

GA - Prineville, Oregon

"I've been having carpal tunnel pains in my hands and wrists for some time. IceWave relieves swelling and pain instantly, with no side effects. I put the white patch on my sternum and the tan patch on which ever hand is giving me a lot of pain. I usually notice the difference in just a few minutes."

DC - Atlanta, Georgia

"I ride my bike to work and found the patches relieved my right shoulder pain completely. My energy level stayed high, I have found I can cope with life emotionally and physically."

DS - Portland, Oregon



LIFEWAVE™
SOFTWARE FOR THE HUMAN BODY